sv BRUKSANVISNING

INSTRUCTIONS FOR USE en

VIKTIGT

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IMPORTANT

Please read all instructions before use.
The product must only be used in accordance with these
instructions. Misuse may endanger the patient's safety.

These instructions for use is also available on the website of the
legal manufacturer www.effectummedical.com/imprint/

Manufacturer:
Effectum Medical AG
Kirchgasse 11
CH-4600 Olten
www.effectummedical.com

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Information for use

INTENDED PURPOSE

ReLink Drain is a breakaway connector for nephrostomy drainage catheter use. ReLink Drain separates into two self-sealable parts when an excessive tension is exerted on the nephrostomy drainage tube, thereby avoiding accidental dislodgement of the nephrostomy catheter, spillage of fluids, and patient injury. The parts are reconnectable after disinfection.

Indications for use

Nephrostomy drainage using a luer lock fitted pigtail catheter secured to the patient using film dressing or equivalent.

Contraindications

Do not use ReLink Drain for:

- Biliary applications.

Intended users

ReLink Drain is only intended for use by healthcare professionals.

Patient target group

Patients receiving nephrostomy drainage therapy except for specified contraindications.

Instructions for use

The following information describes how ReLink Drain is to be used. Users must read the instructions for use carefully before using ReLink Drain for the first time to fully understand all features and indications for use.

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User environment

ReLink Drain is intended to be used in hospitals and other healthcare facilities where nephrostomy drainage treatment is performed.

Intended interaction with patient

External use only.

Duration of use

ReLink Drain is a single use device intended to be used for a maximum of 24 hours and then changed to a new one.

Caution and limitations during use

ReLink Drain should be connected to the extension tube and not directly to the patient's nephrostomy catheter.


ReLink Drain should only be connected to luer lock connections.

ReLink Drain should not be transferred from one drainage line to another.


ReLink Drain should not be left on the patient's drainage line after the drainage is finished.

ReLink Drain must not be reconnected more than five times.

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

ReLink Drain shall be used together with an extension tube between 5 - 25 cm in length (*part 4 in system overview above*).

1. Drainage catheter
2. ReLink Drain
3. Drainage bag
4. Extension tube

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

Unpack and handle ReLink Drain using aseptic technique.

Remove the caps on the female connector and male connector.

Connect male connector directly on the drainage bag and female connector to the extension tube (step 1).

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


Align the two parts and push them together with a steady motion. Please note that when the parts are close to being fully connected, a final hard push is required to activate the membranes and open the flow path.

When ReLink Drain is properly connected, the flow path is open, and fluid can pass through. If the fluid does not flow as expected, it could be because ReLink Drain is not properly connected.

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
If ReLink Drain is accidentally disconnected (step 3a) while in use it can, after disinfection (step 4), be reconnected again. If required, one half of ReLink Drain can be disinfected and the other one changed to a new one (step 3b),

If ReLink Drain has already been reconnected five times and is again disconnected, **the entire ReLink Drain must be replaced** with a new one.

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
Check all surfaces of both ReLink Drain parts. If any visible contaminants are present, replace both parts with a new primed ReLink Drain according to the steps above (step 1-2). If no contaminants are visible, the parts can be disinfected.

Disinfect ReLink Drain membranous parts by continuously rubbing the membrane surfaces for 15 seconds on each membrane with 70% isopropyl alcohol (IPA) wipes (step 4a). Create a pointed tip with the wipe to reach the membranes.

Allow a minimum of 30 seconds for the disinfection agent to dry (step 4b) before reconnecting the ReLink Drain parts according to step 2 above.

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ReLink Drain is made from Polycarbonate, ABS and Silicone and should not be used with fluids not compatible with these materials.

Flowrate of ReLink Drain is 230ml/min measured for viscosity of water, 1.0016 mPa-s at 20°C.

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Warnings

- Do not reuse ReLink Drain as this may change its mechanical features and may cause device failure.
- Do not connect ReLink Drain directly to the patient's nephrostomy catheter as this may lead to bacterial infections.
- Do not use clamps or forceps to connect, tighten or disconnect ReLink Drain from the drainage line or the extension tube, as this can lead to overtightening and damage to ReLink Drain.
- Do not use the same ReLink Drain for more than 24h as this may lead to leakage.
- ReLink Drain should not be transferred from one drainage line to another as this may lead to cross contamination.
- ReLink Drain should not be left on the patient's drainage line after the drainage is finished as this may cause leakage.
- Never connect any other devices to ReLink Drain's membranous parts as this may lead to leakage.
- Do not use ReLink Drain if the use-by date has passed or if the packaging or device is damaged as the device may not be sterile.
- ReLink Drain must not be reconnected more than five times as this may cause device failure and leakage.
- ReLink Drain must not be used with non-luer lock nephrostomy drainage accessories as this may affect the treatment.
- ReLink Drain must not be used for biliary catheters as this may affect the flow in the device.
- ReLink Drain must not be reconnected if not disinfected as this may lead to bacterial infections.

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IMPORTANT

Dispose of used ReLink Drain in a manner consistent with established biohazard procedures in accordance with local regulations.

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All serious incidents that occur in relation to the use of the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

Kontakt / Contact information

Effectum Medical AG | Kirchgasse 11 | CH-4600 Olten | info@effectummedical.com | www.effectummedical.com
AR / EC-REP: Interlinked AB | feedback@interlinked.care